


4 x C-Bored for M3 or #6 SHCAP screws  
 4 x M3x0.5  
 on a 70mm x 60mm rectangle  
 12 x 4-40 UNC



Description:

Date: 9/12/23

Dimensions are in inches [mm]

Nano-BioS series

DRAWING IS TO BE USED AS A GUIDE ONLY. DIMENSIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE. THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MAD CITY LABS. ANY REPRODUCTION IN PART OR AS A WHOLE IS PROHIBITED.